

Specifications



REFRACTIVE POWER MEASUREMENT

| | |
|------------------------|---|
| Spherical power | -30.00 D to +25.00 D (at VD = 12.00 mm) |
| Cylindrical power | 0.00 D to ± 12.50 D (at VD = 12.00 mm) |
| Astigmatic axis | 0° to 180° |
| Minimum pupil diameter | 2.0 mm |

KERATOMETRY MEASUREMENT

| | |
|--------------------------|---------------------|
| Corneal curvature radius | 5.00 mm to 13.00 mm |
| Corneal astigmatic axis | 0° to 180° |

INTRAOCULAR PRESSURE MEASUREMENT

| | |
|-------------------|---------------------------------------|
| Measurement range | 1 mm Hg to 60 mm Hg (1 hPa to 80 hPa) |
|-------------------|---------------------------------------|

PACHYMETRY MEASUREMENT

| | |
|-------------------|------------------|
| Measurement range | 300 µm to 800 µm |
|-------------------|------------------|

TOPOGRAPHY MEASUREMENT

| | |
|--------------------------|---------------------|
| Corneal curvature radius | 5.50 mm to 10.00 mm |
| Corneal astigmatic axis | 0° to 180° |

AUXILIARY FUNCTION

| | |
|-------------------------------------|---|
| Interpupillary distance | 20 mm to 85 mm |
| Corneal diameter and pupil diameter | 1.00 mm to 14.00 mm |
| Dry-eye application | Blinking frequency, tear meniscus height, hyperemia, Meibomian glands |

DATA MANAGEMENT

| | |
|-------------------|--|
| Internal database | Integrated SD card |
| Printer | Integrated thermal printer |
| Data output type | 3x USB 2.0, 1x Ethernet, 1x SD card slot, 1x WLAN (not available in all countries) |
| Export format | DCM, XML, CSV, JPG, PDF |

DIMENSIONS AND ELECTRICAL REQUIREMENTS

| | |
|-------------------|--------------------|
| Dimensions WDH | 312 x 491 x 450 mm |
| Weight | approx. 23 kg |
| Voltage | 100 VAC to 240 VAC |
| Frequency | 50/60 Hz |
| Power consumption | 110 VA |
| Laser class | Class 1 |

TOMEY EUROPE
TOMEY GMBH
Wiesbadener Strasse 21
90427 Nuremberg | Germany
+49 911 938 546 2 - 0
+49 911 938 546 2 - 20
info@tomey.de

tomey.de

Follow TOMEY



TOMEY GmbH is the European headquarters of TOMEY Corporation, 2-11-33 Noritakeshimachi Nishi-Ku, Nagoya, 451-0051, Japan



2023/05 - subject to change without notice

Simplicity meets variety

The versatile multifunctional unit for eye examinations can:

- + Determine the objective refraction
- + Visualise the corneal shape
- + Assess the IOP and C.IOP as part of glaucoma screening
- + Measure the central corneal thickness
- + Observe dry-eye

MR-6000

Multifunction Unit



You + eye.
We care.

MR-6000 Multifunction Unit

The MR-6000 delivers a smart combination of five different eye examinations and a Dry Eye observation app. Along with the advantage of automatic alignment, this means that the MR-6000 speeds up your workflow and makes it more efficient.



"IT'S APPEALING TO WORK WITH HIGH-END TECHNOLOGY DEVICES AND LEARN ABOUT THEIR BENEFITS FOR MEDICAL TREATMENT."

Katharina Koriski

TECHNICAL DEPARTMENT
ENGINEER



Dry-eye application
Blinking frequency, tear meniscus height, hyperemia, and Meibomian gland observation all help assess dry eye.



Topography
16 Mire rings examine an area 8 mm in diameter. A number of topography maps, including Fourier analysis, provide a wide range of options for visualising corneal shape.

Wide range alignment by image recognition
Refined auto alignment and auto shot support the operator during the examination.

Measurement cone interchange in just 4 sec
Electrical cone interchange enables rapid switching between the measurement modes.

Pupil and corneal diameter
The diameters of pupil and cornea (WTW) are measured with each topography. Manual correction allows the user to adjust the measurement points.

Five functions in one system
This multifunction unit includes objective refraction, keratometry, topography, tonometry, and pachymetry. In addition, the Dry Eye observation app completes the powerful 5 + 1 system.

Quick refraction mode
The quick refraction mode supports refraction values to be obtained within seconds despite fixation loss: for example, for patients with nystagmus, children, or uncooperative patients.

Corrected IOP
A new generation of air-flow technology delivers a patient-friendly and gentle air puff. The automated IOP correction mode utilizes pachymetry values to correct the IOP.

